

Cleaning and Disinfecting the GlideScope AVL Systems

Caution. Risk of permanent equipment damage. Do not expose the GlideScope AVL video baton or reusable video laryngoscope to temperatures above 60°C.

AVL Single Use

The GlideScope AVL Single Use is a combination of a nonsterile, reusable video baton and a sterile, single-use Stat. Low Level Disinfection is recommended for the AVL video baton after every patient use. High Level Disinfection is required for the video baton when it is visibly soiled.


AVL Reusable

The AVL laryngoscope is a nonsterile, reusable device. It is required that the AVL Reusable laryngoscope be cleaned and disinfected after every patient use using a High Level Disinfectant method.


1. Disconnect the video cable from the monitor.
2. Insert the attached protective cap into the video cable port to protect the electronic connector as shown.
3. Wash the AVL video baton or reusable blade manually using a hospital grade equipment detergent or an enzymatic debridement agent/detergent to remove all foreign material from the surface of the device.
4. The AVL video baton or reusable blade can now be disinfected.



Video cable disconnected from monitor



Cap inserts for cleaning



Reusable blade

Video baton

Cleaning and Disinfecting the Monitor

Clean the exterior of the monitor by wiping with IPA (70% isopropyl alcohol), bleach (100ppm), or a mild detergent and water.



Monitor front

Cleaning and Disinfecting the Cradle

Clean the cradle by wiping with a standard hospital-grade surface cleaning product.



Cradle

Risk Assessment

Device	Sterile	Use	Spaulding's/ CDC Classification	Disinfection			Sterilization
				Low Level	Intermediate Level	High Level	
GVL® Stat	Sterile	Single Use	Semicritical				
Video Baton*	Non-sterile	Reusable	Noncritical	X			
AVL Reusable	Non-sterile	Reusable	Semicritical			X	
GlideRite® Rigid Stylet	Non-sterile	Reusable	Semicritical			X	
Monitor#	Non-sterile	Reusable	Noncritical	X			
Cradle#	Non-sterile	Reusable	Noncritical	X			

IT IS UNDERSTOOD THAT ALL ITEMS IN THIS CHART WILL BE USED AS INTENDED.

Key:

* The AVL video baton is a nonsterile reusable device, which is protected from contact with mucous membranes and nonintact skin by the Stat (sterile, single-use) when used as intended. Low Level Disinfection is recommended for the AVL video baton after every patient use. High Level Disinfection is required for the video baton when it is visibly soiled.

Disinfect monitor and cradle when they are visibly soiled and on a regular basis as per a schedule established by the medical care facility or provider.

- Shaded areas show levels of disinfection that are not required/not compatible with device materials.
- Checked boxes (x) show minimum requirement.
- Unshaded areas show permissible levels of disinfection based on compatibility with device materials.

Warning: Disinfectants and cleaning methods listed are recommended based on compatibility with device materials. Refer to the disinfectant label instructions for guidance on disinfection efficacy and appropriate clinical uses.

Caution: Meticulous cleaning must precede any disinfection process, to ensure all foreign matter is removed from the surface of the device. This allows the active ingredients of the chosen process to reach all the surfaces of the device.

Note: When using any high-level disinfectants, read and comply with product use instructions in all applications.

Availability of disinfection products varies by country, and we are unable to test products in every market. Please use the list of recommended disinfectants in this manual to compare with products available locally.