

Advanced video laryngoscopy for small patients has never been more clear.

The GlideScope AVL defines advanced video laryngoscopy—airway views in DVD clarity, with real-time recording to capture the details of your difficult airway cases.

With a digital color monitor, the GlideScope AVL gives you four single-use size options for children from preterm to 28 kg—all provide a consistently clear view, enabling swift intubation.

There's also a newly designed, reusable GVL® 2 that features a slimmer profile for small patients.

AVL=ADVANCED AIRWAY VIEWS

- ▶ Digital color monitor with DVD-quality resolution
- ▶ Reusable video baton or GVL 2 with high-resolution camera
- ▶ Advanced resolution output for OR theater viewing
- ▶ Anti-fogging mechanism resists lens clouding and secretions

AVL=ADVANCED USER FEATURES

- ▶ Integrated, real-time recording helps confirm tube placement and facilitate teaching
- ▶ Onboard *GlideScope 4-Step Technique* tutorial makes AVL easy to use, learn, and teach
- ▶ Ergonomic design and rugged construction
- ▶ 4 single-use Stat sizes, in sterile packaging
- ▶ 1 reusable GVL blade for infant to small child



Advanced Video Laryngoscopy

Real-time recording. DVD clarity.



Typical airway view.






Applications


- First-use intubations, replacing direct laryngoscopy
- Normal or restricted oropharyngeal views
- Cormack-Lehane grades I-IV laryngeal views
- Trauma airways—excellent when dealing with blood and secretions in the airway
- Nasal tracheal intubation
- Video-guided tube exchange in the NICU



System Components

- Digital, DVD-quality color monitor
- Select from reusable video baton or newly designed, reusable GVL® 2
- Single-use sterile Stat, choice of 4 sizes (box of 10)
- Choice of mobile stand or IV pole mount
- User’s Manual and Quick Reference Cards
- 1-year warranty

GlideScope Products and Systems—Designed for 1st Pass Success

SINGLE-USE SYSTEM				
AVL VIDEO BATON	GVL 0	GVL 1	GVL 2	GVL 2.5
Video Baton 1-2	< 1.5 kg	1.5 - 3.6 kg	1.8 - 10 kg	10 kg - 28 kg
 <p>Length of flexible baton: 66 mm Height at camera: 6 mm Width at camera: 7 mm Cable length: 204 cm Weight: 180 g</p>	 <p>Length (blade tip to handle): 36 mm Height at camera: 9 mm Width at camera: 9 mm Max blade width in front of camera: 9 mm</p>	 <p>Length (blade tip to handle): 44 mm Height at camera: 9 mm Width at camera: 10 mm Max blade width in front of camera: 13 mm</p>	 <p>Length (blade tip to handle): 56 mm Height at camera: 9 mm Width at camera: 11 mm Max blade width in front of camera: 16 mm</p>	 <p>Length (blade tip to handle): 63 mm Height at camera: 9 mm Width at camera: 13 mm Max blade width in front of camera: 20 mm</p>

REUSABLE SYSTEM
GVL 2
4 - 20 kg
 <p>Length (blade tip to handle): 51 mm Height at camera: 11 mm Width at camera: 14 mm</p>

STANDS	
MOBILE STAND	IV POLE MOUNT
 <p>A stand-alone configuration for use in a variety of hospital settings.</p> <p>Power source: Battery and/or AC power Monitor weight: 1.0 kg Mobile stand weight: 3.2 kg</p>	 <p>Easily integrates into all standard hospital environments.</p> <p>Monitor weight: 1.0 kg IV pole mount weight: .9 kg</p>

Note: System configurations vary.

GlideScope, the GlideScope symbol, GVL, Verathon, and the Verathon Torch symbol are trademarks of Verathon Inc. © 2010, 2011 Verathon Inc.

GlideScope video laryngoscope systems are CE marked in accordance with the Medical Device Directive, and the Verathon Inc. quality system is Quality System Certified to ISO 13485:2003 standards.

0900-2699-04-86



Corporate Headquarters

20001 North Creek Parkway / Bothell, WA 98011
425.867.1348 / 800.331.2313 / 425.883.2896 Fax
verathon.com