

For advanced video laryngoscopy, the choice has never been more clear.

The GlideScope AVL Single Use features a digital color monitor and digital camera, providing a consistently clear, DVD-quality view, enabling swift intubation.

The AVL gives you single-use options designed for a wide range of patient weights and sizes—from preterm to morbidly obese—ideal for demanding OR, ED, and ICU settings.

Add real-time recording and an onboard video tutorial, and these versatile products are a clear choice in advanced video laryngoscopy.



#### AVL=ADVANCED AIRWAY VIEWS

- ▶ Digital color monitor with DVD-quality resolution
- ▶ Reusable video baton with digital-resolution camera
- ▶ Advanced resolution output for OR theater viewing
- ▶ Anti-fogging mechanism resists lens clouding and secretions

#### AVL=ADVANCED USER FEATURES

- ▶ Integrated, real-time recording helps confirm tube placement and facilitate teaching
- ▶ Onboard *4-Step Technique* tutorial makes AVL easy to use, learn, and teach
- ▶ Ergonomic design and rugged construction for multiple clinical applications
- ▶ 6 sterile, single-use Stats, designed for preterm to morbidly obese patients



## Advanced Video Laryngoscopy: Real-time recording. DVD clarity.



### Applications

- First-use intubations, replacing direct laryngoscopy
- Normal or restricted oropharyngeal views
- Cormack-Lehane Grades I–IV laryngeal views
- Trauma airways—excellent when dealing with blood and secretions in the airway
- Nasal tracheal intubation
- Video-guided tube exchange in the ICU

### System Components

- Digital, DVD-quality color monitor
- Reusable video baton with digital camera
- Single-use Stat, choice of 6 sizes (box of 10)
- GlideRite® Rigid Stylet
- Protective video baton cradle
- Choice of mobile stand or IV pole mount
- User’s Manual and Quick Reference Cards
- Instructional In-Service video

## GlideScope Products and Systems: Designed for 1<sup>st</sup> Pass Success

		SINGLE USE STATS					
GVL® Size		GVL 0	GVL 1	GVL 2	GVL 2.5	GVL 3	GVL 4
Patient Weight		<1.5 kg	1.5–3.6 kg	1.8–10 kg	10–28 kg	10 kg–Adult	40 kg–Morbidly Obese
AVL	Video Baton 1-2						
	Video Baton 3-4						
GlideRite® Rigid Stylet*							

AVL System	AVL Video Batons	GVL 0, 1, 2, and 2.5	GVL 3 and 4
	<p><b>Video Baton 1-2</b>                      Length of flexible baton: 66mm                      Height at camera: 6mm                      Width at camera: 7mm                      Video cable length: 204cm                      Weight: 180g</p> <p><b>Video Baton 3-4</b>                      Length of flexible baton: 105mm                      Height at camera: 11mm                      Width at camera: 11mm                      Video cable length: 203cm                      Weight: 238g</p>	<p>Length (blade tip to handle): 36/44/56/63mm                      Height at camera: 9/9/9/9mm                      Width at camera: 9/10/11/13mm                      Max blade width in front of camera: 9/13/16/20mm</p>	<p>Length (blade tip to handle): 78mm/92mm                      Height at camera: 16mm/15mm                      Width at camera: 16mm/20mm                      Max blade width in front of camera: 21mm/27mm</p>

Note: System configurations vary.  
 \*For use with endotracheal tubes 6.0 mm and larger.

AVL PRODUCT CONFIGURATION	PRODUCT OPTIONS
<p><b>On Mobile Stand</b></p> <p>A stand-alone configuration for use in a variety of hospital settings</p> <p>Power source: battery and/or AC power</p> <p><b>Also available in IV Pole Mount</b></p>	<p><b>GlideScope Direct</b></p> <p>The GlideScope Direct intubation trainer is designed to facilitate instruction of classic direct laryngoscopy (DL).</p> <p><b>Available in Single-use and Reusable</b></p>

GlideScope, the GlideScope symbol, GVL, GlideRite, Verathon, and the Verathon Torch symbol are trademarks of Verathon Inc. © 2010 Verathon Inc.

GlideScope® video laryngoscope systems are CE marked in accordance with the Medical Device Directive, and the Verathon Inc. quality system is Quality System Certified to ISO 13485:2003 standards.

0900-2695-04-86



**Corporate Headquarters**  
 20001 North Creek Parkway / Bothell, WA 98011  
 425.867.1348 / 800.331.2313 / 425.883.2896 Fax  
 verathon.com